



# KFU

## جامعة الملك فيصل

KING FAISAL UNIVERSITY

جامعة ووطن.. نماء.. واستدامة..



---

# INSTRUMENT STANDARD OPERATING PROCEDURE MANUAL

---

College of Medicine



SAFETY AND LABORATORY COOMMITTEE, C.O.M, K.F.U.

**Prepared by**

Document Number	Name	Signature	Date
CM 001 CM/BCH-041	GENESYS™ 10 Series Spectrophotometers		
Revision Number	Approved by		
R0 Write the revision number	Dean:  Department:		

**Document History**

Document Section	Details of Amendments	Date	Modified by (Initials)
SOP	First Draft on SOP for the operation of ( GENESYS™ 10 Series Spectrophotometers)		

**1. OBJECTIVE**

- ☐ The document describes the operation of ( GENESYST™ 10 Series Spectrophotometers)

## 2. SCOPE

The spectrophotometer, shown in Figure 1, is an

easy-to-use, UV-Visible spectrophotometer  
designed for quantitative measurements in:

- ☐ •Industrial quality control and research  
Laboratories, including food and beverage,  
Chemical and pharmaceutical applications
  - ☐ •Academic research and teaching laboratories
  - ☐ •Environmental laboratories, including water  
and wastewater applications
  - ☐ •Life science laboratories, for working on  
Proteins and nucleic acid research
- Your spectrophotometer offers:
- ☐ •Unique optical system that ensures accuracy  
and precision
  - ☐ •Automatic wavelength calibration on power-up  
to ensure that the instrument is operating  
properly
  - ☐ •Built-in applications software for concentration,  
standard curve, absorbance ratio, absorbance  
difference, survey scanning, kinetics,  
multi-wavelength analysis and three-point net  
absorbance
  - ☐ •Built-in procedures that allow you to validate  
the performance of your instrument.

## 3. RESPONSIBILITIES

- ☐ It is the responsibility of designated personnel in the lab to train staff and students  
on this procedure and to ensure adherence to this procedure under supervision.
- ☐ It is the responsibility of designated personnel (staff or Student) to follow the  
instructions of this procedure under supervision.

## 4. REFERENCES

- ☐ There is an Operating Manual provided by the Manufacturer.

## 5. DEFINITIONS

- ☐ Mention the full form of any abbreviations used.
- ☐ Mention the standard definition of any term used.

## 6. SAFETY PRECAUTIONS

- ☐ **Line voltages**  
100 - 240 VAC  $\pm 10\%$   
50 - 60 Hz  
80 VA max
- ☐ **Operating environment**  
The instrument meets the specifications on the previous page under the following conditions  
After a 30-minute warm-up period.
- ☐ **Ambient temperature:** 5° to 35° C  
(41 °F to 95° F)
- ☐ **Relative humidity:** 20% - 80% RH
- ☐ **Storage environment**  
-20 °C to 70° C (-4° F to 158° F) Relative humidity not to exceed 0.040 pounds moisture per pound of dry air. Allow instrument to adjust to room temperature for 24 hours after taking it out of storage.
- ☐ **Temperature should be maintained at  $\pm 4^\circ$  F.**  
Relative humidity should be maintained to  $\pm 5\%$ ..

## 7. PROCEDURE FOR OPERATING ( GENESYS™ 10 Series Spectrophotometers)

### 7.1. Turning on the instrument:

Press the back button, and the device will operate

- ☐ **Software Operation.**

### 7.2. Steps of the Procedure

Wait for the data to appear on the LCD

Select the job you want to do:

For example: (measuring absorption).

1 -Zero the device by pressing the plank gauge, lift your finger and wait for the zero to appear on the LCD.

2 -Enter the wavelength at which you want to measure the absorption.

3 -Lift the cap and place the sample that you want to measure in the place provided for it on the right side, then close the cover.

4- Wait for the result to appear on the LCD.

### **7.3. Turning off the Instrument**

Press the back button to turn off the device.

### **7.4. Warning**

**Write any the warning in BOLD and RED FONT.**

---