



KFU

جامعة الملك فيصل

KING FAISAL UNIVERSITY

جامعة ووطن.. نماء.. واستدامة..

INSTRUMENT STANDARD OPERATING PROCEDURE MANUAL

College of Medicine



SAFETY AND LABORATORY COOMMITTEE, C.O.M, K.F.U.

Prepared by

Document Number	Name	Signature	Date
CM/CTC-001	Mohammad Ahmad Al-oudah		
Revision Number	Approved by		
R0	Dean:		
Write the revision number	Department:		

Document History

Document Section	Details of Amendments	Date	Modified by (Initials)
SOP	Automated Delivery NOELLE S2200		

1. OBJECTIVE

Training students to detect newborns and all possible cases.

2. SCOPE

Training medical students on natural and unnatural birth operations.

3. RESPONSIBILITIES

- ☐ It is the responsibility of designated personnel in the lab to train staff and students on this procedure and to ensure adherence to this procedure under supervision.
- ☐ It is the responsibility of designated personnel (staff or Student) to follow the instructions of this procedure under supervision.

4. REFERENCES

Device Catalog

5. DEFINITIONS

- ☐ Mention the full form of any abbreviations used.
- ☐ Mention the standard definition of any term used.

6. SAFETY PRECAUTIONS

Electric device running all the time

7. PROCEDURE FOR OPERATING write Instrument Name)

7.1. Turning on the instrument:

Obstetrics and Gynecology

- 1- Mother doll (large)
- 2- Fetal doll (electronic).
- 3- Baby doll (regular)
- 4- Tablet (controller)
- 5- A computer (run the program and display the vital signs of the mother).
- 6- Screen display device (to display the vital elements of the fetus).

7.2. Software Operation.

The program is located on the desktop, the controller and the monitor with the name Gaumard

7.3. Steps of the Procedure

Charge mother doll and baby doll

Turn on the tablet and the computer and choose the VICTORIA program from the desktop in both the device that is running the computer, to connect it to the tablet with Wi Fi

We open the mother's belly through buttons located in the middle of her waist and on left

We put a sufficient amount of oil for the doll on the internal parts of the mother that are in contact with the fetus and the fetus itself, as well as the placenta (to .preserve the mother's organs and nostalgia from rupture during childbirth)

The fetus is placed on a special platform to stabilize it and a stent is installed for it in its back. Then the fetus is placed in. The stent is attached to the thrusting machine in the fetus, and the placenta is attached to the inner wall of the mother's .abdomen

After the birth process, we open the mother's abdomen through the buttons in the .position of the fetus, and so on in each operation

7.4. Turning off the Instrument

After completing the practical lesson, there is an (X) at the top of the screen on the right and we choose to save the session or not.

7.5. Warning

Electric device running all the time

.



Safety and Laboratory Committee, CO