



KFU

جامعة الملك فيصل

KING FAISAL UNIVERSITY

جامعة ووطن.. نماء.. واستدامة..



INSTRUMENT STANDARD OPERATING PROCEDURE MANUAL

College of Medicine



SAFETY AND LABORATORY COOMMITTEE, C.O.M, K.F.U.

Prepared by

Document Number	Name	Signature	Date
CM 001	MR. HUSAIN AL NWAISER MR. ALI AL THANI		21.4.2021
Revision Number	Approved by		
R0	Dean: Department:		

Document History

Document Section	Details of Amendments	Date	Modified by (Initials)
SOP	First Draft on SOP for the operation of Water Purification Systems(Milli-Q)	21.4.2021	

1. OBJECTIVE

- ☐ The document describes the operation of **Water Purification Systems(Milli-Q)**.
- ☐ **Uses regular tap water as feed, and produces pure and ultrapure water delivered by independent POD .**
- ☐ **Allows control over both ionic and organic contaminants that can impact research results.**
- ☐ Allows installation of a comprehensive water purification chain in a minimal amount of space.
- ☐ Eliminates the constraints of centralized water stations and their distribution loops (complexity, heavy maintenance, prolonged downtime, contamination)
- ☐ To obtain Milli-Q water quality, purified water passes through a long-life dual wavelength UV lamp, which ensures organic molecule degradation by photooxidation.– Ultrapure water free of particulates, bacteria, pyrogens, nucleases, proteases, VOCs, endocrine disruptors and organics .
- ☐ Feed water quality Elix, RiOs, distilled or DI water with conductivity < 100 $\mu\text{S}/\text{cm}$ and TOC < 50 ppb Feed water connection 1/2 in Gaz M Feed water pressure.
- ☐ Electrical power supply frequency 50 – 60 Hz \pm 10 % Feed water connection 1/2 in. Gaz F Data connection Ethernet (RJ45)* 0 – 0.3 bar Feed water temperature 5 – 35 °C
- ☐ The Main Display is used to navigate the Milli-Q System software.

2. SCOPE

3. RESPONSIBILITIES

- ☐ It is the responsibility of designated personnel in Research lab to train staff and students on this procedure and to ensure adherence to this procedure under supervision.
- ☐ It is the responsibility of designated personnel (staff or Student) to follow the instructions of this procedure under supervision.
- ☐ Routine in R & D laboratories of the pharmaceutical and biotechnical industry.

4. REFERENCES

- ☐ **Water Purification Systems(Milli-Q)** operating instructions manual

5. DEFINITIONS

- | | |
|------------------------------|---------------------------------------|
| 1-Tap water feed | 7-Feed water conductivity cell |
| 2-Strainer | 8-RO cartridge with sanitization port |
| 3-Inlet solenoid valve | 9-RO reject recycling |
| 4-Progard® pretreatment pack | 10-Drain |
| 5- Pressure regulator | 11-RO reject solenoid valve |
| 6- Booster pump | 12-Permeate conductivity cell |

13-3-way solenoid valve	22-E-POD® dispenser and final polisher
14-Elix® electrodeionization module	23-Photooxidation UV Lamp
15-Elix® resistivity cell	24-Quantum® cartridge
16-Thermistor	25-Milli-Q® product water resistivity cell
17-Check valve	26-Thermistor
18-Bactericidal UV Lamp	27-Milli-Q® water TOC monitor
19-PE Reservoir, ASM and Vent filter	28-Q-POD® dispenser
20-Delivery pump	29-Final polisher Application-Pak
21-Flow meter	30-Recirculation solenoid & check valve

6. SAFETY PRECAUTIONS

- ☐ Do not remove the covers of the Milli-Q System at any time. Electrical and mechanical components inside the Milli-Q System could pose a hazard.
- ☐ A qualified Millipore Service Representative should perform any work that needs to be done while the Milli-Q System is opened.
- ☐ Make sure the PS/2 cable and the Termination Plug are well connected to the Q-POD Unit.

7. PROCEDURE FOR OPERATING THE Water Purification Systems(Milli-Q)

7.1. Turning on the instrument:

- ☐ Switch on the instrument (mains switch on).

7.2. Software Operation.

7.3. Steps of the Procedure

- ☐ After setup and installation .
- ☐ Switch start up /On

7.4. Turning off the Instrument

- ☐ Switch the unit off and disconnect from the mains

7.5. Warning

Write any the warning in BOLD and RED FONT.
