



KFU

جامعة الملك فيصل

KING FAISAL UNIVERSITY

جامعة ووطن.. نماء.. واستدامة..



INSTRUMENT STANDARD OPERATING PROCEDURE MANUAL

College of Medicine



SAFETY AND LABORATORY COOMMITTEE, C.O.M, K.F.U.

Prepared by

Document Number	Name	Signature	Date
CM/EMU-001	husain musa al nwaiser ali althani		
Revision Number	Approved by		
R0	Dean:		
Write the revision number	Department:		

Document History

Document Section	Details of Amendments	Date	Modified by (Initials)
SOP	First Draft on SOP for the operation of Transmission Electron Microscopy (TEM- JEOL(JEM-1011)		

1. OBJECTIVE

The document describes the operation of Transmission **Electron Microscopy (TEM- JEOL (JEM-1011))**.

- ☐ Used to study the internal structures of the cells.

2. SCOPE

- ☐ Full automated software processor system with integrated DIGITAL CAMERA.
- ☐ Accelerating voltage 40 to 100Kv.
- ☐ Magnification mode SAP 10B :800- 600,000 TIMES
- ☐ For biological samples (virus, tissue...etc and Geological specimen)

3. RESPONSIBILITIES

- ☐ It is the responsibility of designated personnel in the lab to train staff and students on this procedure and to ensure adherence to this procedure under supervision.
- ☐ It is the responsibility of designated personnel (staff or Student) to follow the instructions of this procedure under supervision.

4. REFERENCES

- ☐ FOR Transmission Electron Microscopy instructions manual (TEM- JEOL(JEM-1011)).

5. DEFINITIONS

- ☐ Transmission Electron Microscopy (TEM- JEOL(JEM-1011)) IS A Closed fully automated controlled programming software..
- ☐ Control Panel – R (FOR MAGNIFICATION, FOCUSING, FILAMENT AND CONTROL HIGH TENSION).
- ☐ Control Panel - L (Adjust the contrast and brightness.).
- ☐ Control Panel – CL (Switch on and open the Right control panel).
- ☐ Control PF (to specify the voltage maximum 100 KV)
- ☐ Column.
- ☐ Four Electron magnetic Lenses. (condenser, objective, intermediate and projector).
- ☐ VACUUM TUBE.
- ☐ Chiller.

- ☐ DP(Air Power Vacuum).
- ☐ Specimen holder.(2 specimen each)
- ☐ Specimen Section
- ☐ Electron beam(A node).
- ☐ Electron Gun.
- ☐ Fluorescent Screen.
- ☐ Digital camera(0.5 nm).
- ☐ Air compressor.
- ☐ Screen
- ☐ Printer

6. SAFETY PRECAUTIONS

- ☐ Do not remove the grounding wire ,it may cause electric shock.
- ☐ Do not dismount or reassemble thhe microscope column.
- ☐ Any abnormalities occur immediately shut down and contact the JEOL ENGIENEER .
- ☐ Do not touch the Wehnrlt to avoid the electric shock.
- ☐ Do not put your hand between the upper and lower electron gun chambers.
- ☐ Before starting be sure the gun lift and lens power set as ON.

6.1. PROCEDURE FOR OPERATING Transmission Electron Microscopy (TEM-JEOL(JEM-1011)) IS ATurning on the instrument:

6.2. Software Operation.

6.3. Steps of the Procedure

- ☐ **Turning on the regulation power supply (Leave it overnight).**
- ☐ Turn On the control panels (L-R,CL,PF).
- ☐ Open the compressed air pressure(0.28-0.47)MPa.
- ☐ Open the main valve of cooling water supply .
- ☐ Check the personal file for the operating condition.
- ☐ Place the specimen holder onto the specimen holder stand.
- ☐ Loading the films into a dispensing magazine. And place the magazine be into of camera chamber
- ☐ Set the voltage value for the electron beam – for biologic samples is 80Kv (maximum voltage is 100 Kv
- ☐ Align the condenser lenses.

- ☐ Inserting the CONDENSER LENSE into the beam path.
- ☐ Condenser lens astigmatism correction
- ☐ Obtain the filament image.
- ☐ Adjust the specimen by using goniometer.

6.4. Turning off the Instrument

END OF PROGRAM

- ☐ **To turn off the instrument.**
- ☐ Unload the samples
- ☐ Turn off the software.
- ☐ Turn off the control panel..
- ☐ Turn off the electric main power source.
- ☐ Wait for 15 minutes
- ☐ Cut the power from the control panel.

6.5. Warning

Write any the warning in BOLD and RED FONT.

