



KFU

جامعة الملك فيصل

KING FAISAL UNIVERSITY

جامعة ووطن.. نماء.. واستدامة..



INSTRUMENT STANDARD OPERATING PROCEDURE MANUAL

College of Medicine



SAFETY AND LABORATORY COOMMITTEE, C.O.M, K.F.U.

Prepared by

Document Number	Name	Signature	Date
CM/PHY-063	Mr.Ibrahim Al-Saqer Mr.AbdulRazaq Alwebari		11-02-2022
Revision Number	Approved by		
R1	Dr.Tarek BENAMEUR Department:Physiology		

Document History

Document Section	Details of Amendments	Date	Modified by (Initials)
SOP	First Draft on SOP for the operation of Ultracentrifuge Optima MAX-XP		

1.OBJECTIVE:

- ☐ The document describes the operation of Ultracentrifuge Optima MAX-XP

2.Description:

The Optima MAX-XP microprocessor-controlled tabletop ultracentrifuge generates high centrifugal forces for a variety of applications. The ultracentrifuge design features a variable frequency induction drive, thermoelectric temperature control system, self-purging vacuum system, rotor overspeed identification system, user login feature, program memory that contains multiple five-step programs, and a choice of acceleration and deceleration rates. Manual and programmed operations are available from the integrated touchscreen interface.

- In manual operation, you enter the individual run parameters before beginning each run.
- In programmed operation, you can duplicate runs quickly and accurately by selecting previously entered programs and running them again.

3.RESPONSIBILITIES:

- ☐ It is the responsibility of designated personnel in the lab to train new staff and students on this procedure and to ensure adherence to this procedure under supervision.
- ☐ It is the responsibility of designated personnel (staff or Student) to follow the instructions of this procedure under supervision.
- ☐ The head of the physiology department must resolve any problem with the process and difficulties in using this SOP.

4.REFERENCES:

- ☐ Ultracentrifuge Optima MAX-XP Instruction Manual, BECKMAN COULTER

5.SAFETY PRECAUTIONS:

- ☐ Wear gloves, a lab coat and safety glasses for protection.

6.PROCEDURE FOR OPERATING of Ultracentrifuge Optima MAX-XP

6.1.Turning on the instrument:

Turn the power on for the device.

6.2.Creating a New Program:

- 1- In the Main screen, press the MENU button.
- 2- Select Program.
- 3- Press the NEW PROGRAM button
- 4- Select a rotor first
- 5- Select the desired rotor, and press the OK button to accept.
- 6- Press the NEW STEP button.
- 7- Enter the speed on the keypad display, and press the OK button to accept.
- 8- Enter the run time on the keypad display, and press the OK button to accept.
- 9- Enter the temperature on the keypad display, and press the OK button to accept.
- 10- Repeat Steps 6–9 to enter the parameters for up to five steps for a complete run.
- 11- To change the acceleration/deceleration rates, press either the ACCEL or DECEL button.
- 12- Select the desired rates by touching the corresponding numbers.

- 13- Press the OK button to accept.
- 14- Press the OK button to accept.
- 15- Type the name of the program using the keypad just as you would a keyboard.
- 16- press the OK button to accept.
- 17- Repeat this procedure to add and save additional programs.

6.3.Starting a Programmed Run:

- 1- In the Main screen, press the MENU button.
- 2- Select Program
- 3- Select the name of the program you want to run, and press the OK button.
- 4- Press the START button.
- 5- To terminate a run for any reason, press the STOP button
- 6- To exit from program mode:
 - a. Press the MENU button in the Main screen.
 - b. In the Select Program screen, select No Program Selected.
 - c. Press OK.
 - The Main screen returns to view so that you can enter parameters manually.
 - d. Alternatively, you can try changing a parameter such as speed, time, or temperature.
 - A message appears asking you to confirm that you want to exit from program mode.

6.4.Turning off the Instrument:

- ☐ Turn the power off for the device.

6.5.Warning:

All patients specimens should be considered potentially infectious and must be handled with precautions used for human blood
