



KFU

جامعة الملك فيصل

KING FAISAL UNIVERSITY

جامعة ووطن.. نماء.. واستدامة..



INSTRUMENT STANDARD OPERATING PROCEDURE MANUAL

College of Medicine



SAFETY AND LABORATORY COOMMITTEE, C.O.M, K.F.U.

Prepared by

Document Number	Name	Signature	Date
CM/CTC-069	Mohammed Ahmad AL-oudah		
Revision Number	Approved by		
R0	Dean:		
Write the revision number	Department:		

Document History

Document Section	Details of Amendments	Date	Modified by (Initials)
SOP	LAP X		

1. OBJECTIVE

Laparoscopic surgery simulator

2. SCOPE

Training students to acquire basic skills of laparoscopic surgery

3. RESPONSIBILITIES

- ☐ It is the responsibility of designated personnel in the lab to train staff and students on this procedure and to ensure adherence to this procedure under supervision.
- ☐ It is the responsibility of designated personnel (staff or Student) to follow the instructions of this procedure under supervision.

4. REFERENCES

Device Catalog

5. DEFINITIONS

- ☐ Mention the full form of any abbreviations used.
- ☐ Mention the standard definition of any term used.

6. SAFETY PRECAUTIONS

Electric device running all the time

7. PROCEDURE FOR OPERATING write Instrument Name)

7.1. Turning on the instrument:

- 1- A computer (with two Display Screen) the program is installed on the desktop
- 2- Four of the laparoscopic surgery instruments by hand simulator

7.2. Software Operation.

The program is available on the desktop only

7.3. Steps of the Procedure

- 1- Turn on the computer
- 2- Choose the Simulator-of-laparoscopic-surgery program from the desktop
- 3- Enter MODE USER: kfuadmin
Password: medical
- 4- Follow the instructions on the device screen to configure the simulator of laparoscopic surgery tools
- 5- Exit the list that contains all the operations be trained

7.4. Turning off the Instrument

After completing the training go to the main menu and choose SELECT USER

MOD USER: 427 then you are a common no password required.

(X) Mark appears at the top of the program page on the right, and after that, you log out

7.5. Warning

8. **Electric device running all the time**

