



KFU

جامعة الملك فيصل

KING FAISAL UNIVERSITY

جامعة ووطن.. نماء.. واستدامة..



INSTRUMENT STANDARD OPERATING PROCEDURE MANUAL

College of Medicine



SAFETY AND LABORATORY COOMMITTEE, C.O.M, K.F.U.

Prepared by

Document Number	Name	Signature	Date
CM/SIM-030	Rami ALgamdi		
Revision Number	Approved by		
R0	Dean:		
Write the revision number	Department:		

Document History

Document Section	Details of Amendments	Date	Modified by (Initials)
SOP	P.A.C. Complete Camera System		

1. OBJECTIVE

P.A.C. Complete Camera System

2. SCOPE

Training students on how to use knee and shoulder endoscopes and gain training the skill of dealing with them

3. RESPONSIBILITIES

- ☐ It is the responsibility of designated personnel in the lab to train staff and students on this procedure and to ensure adherence to this procedure under supervision.
- ☐ It is the responsibility of designated personnel (staff or Student) to follow the instructions of this procedure under supervision.

4. REFERENCES

Device Catalog

5. DEFINITIONS

- ☐ Mention the full form of any abbreviations used.
- ☐ Mention the standard definition of any term used.

6. SAFETY PRECAUTIONS

Electric device running all the time

7. PROCEDURE FOR OPERATING write Instrument Name)

7.1. Turning on the instrument:

1-Camera

2-Arthroscopy Scope

3-Light Cable

4- Light Source

5-P A C Camera Installation CD

7.2. Software Operation.

Installing the software on a computer (it is on CD)

7.3. Steps of the Procedure

Install the program on a computer and we install the camera with the arthroscopy scope and connect it with the light source by a light cable. Run the program from the desktop in the work of the device

7.4. Turning off the Instrument

Log out of the program page

7.5. Warning

Contains breakable parts of the camera and optical fibers (optical transmission cable)

